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Prefabricated hybrid ceramic blocks, compatible with all dental CAD/CAM systems.

DESCRIPTION

NUMERYS HC, as a pre-cured hybrid ceramic block type, is a dental restorative material used to make inlays, onlays, veneers and crowns. Our blocks are very easy to apply to the appropriate patient teeth's shade and to use on dental CAD/CAM's milling procedures.

COMPOSITION

The monomer matrix is composed of Bis-GMA and TEGDMA. The inorganic fillers are barium aluminosilicate, silicon dioxide and zirconia. The total content of inorganic fillers is ca.80wt% and the particle size of inorganic fillers is between 10nm and 3µm. Additives are initiators, stabilizers, pigments...

INDICATIONS

Fabrication of inlays, onlays, veneers and crowns with radiopaque properties.

RECOMMENDED PROCEDURE

Machining

1. Prepare the tooth according to existing standards
2. Select the appropriate block's shade
3. Select the block's size according to the restoration
4. Select the suitable machining program for hybrid ceramic blocks
5. Insert the mandrel without forcing and mill the block to the proper size using a suitable CAD/CAM system

Bonding

1. Test the machined prosthesis in mouth to, if necessary, adjust it with polishers and diamond burs
2. Before bonding, apply a silane primer (we recommend Silanea)
3. Apply a permanent self-adhesive or dual-cure resin cement (we recommend Total C-RAM) and place the definitive restoration on the tooth/teeth
4. If applicable, light-cure during 2 sec, remove the cement's excess and finish by polishing (we recommend Perfect Polish)
5. If necessary, stain the prosthesis with suitable dyes (we recommend light-cured characterization materials)

STORAGE

1. Avoid exposure to the sunlight
2. Store at room temperature

PRECAUTIONS

1. NUMERYS HC should be used by a dentist or dental professionals
2. Keep away from a reach of non-dental users
3. Use of protective gloves and glasses is recommended
4. Do not use after the expiry date
5. If skin sensitization or allergy (methacrylate for ex.) occurs, stop the use of this product
6. Do not store or expose at elevated temperature or intense light
7. Apply for dental purpose only
8. Be careful not to fall down the product, especially on the mandrel
9. If contamination occurs, do not use the product

WARRANTY

ITENA CLINICAL will replace any product that is proven to be defective. ITENA CLINICAL does not accept liability for any damage or loss, direct or consequential, stemming from the use of or inability to use the product described. Before using it, it is the responsibility of the user to determine the suitability of the product for its intended use and the suitable machining program for hybrid ceramic. The user assumes all risks and liability in connection therewith.

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